

Aluminum L-Lactate

Art. ID TRC-A575730-5G
Unit 5 g
Deliverydetails /not restricted

Description

Category: Miscellaneous Reagents /// Application Notes: Aluminum L-Lactate has been used to inject aluminum into rabbits. Aluminum L-Lactate has been shown to decrease human enamel dissolution in a pH 4 acetate buffer. /// References: Melograna, J., Yokel, R. A.: Res. Commun. Chem. Path., 44, 411 (1984)| Putt, M. S., Kleber, C.J.: J. Dent. Res., 64, 437 (1985)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminum L-lactate	[18917-91-4]				