

Aluminum Isopropoxide

Art. ID TRC-A575930-1KG
Unit 1 kg
Deliverydetails UN3181 II

Description

Category: Miscellaneous Reagents /// Appearance: Solid /// Application Notes: Aluminum isopropoxide is a compound used as an initiator for the polymerization of various compounds, such as lactides. Aluminum isopropoxide is also used as a reducing agent for the reduction of Cinnamaldehyde (C442020). /// References: Dubois, P., et al.: Macromolecules. 24, 2266 (1991)| H?jek, J., et al.: Appl. Catal. A: Gen., 251, 385 (2003)| Lou, X., et al.: Macromol. Rapid Comm., 23, 126 (2002)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminum isopropoxide	[555-31-7]				