

(R,R)-trans-3-Amino-1-indanol

Art. ID TRC-A612120-250MG
Unit 250 mg

Description

Category: Amines, Aromatics /// Application Notes: (R,R)-trans-3-Amino-1-indanol is an isomer of 3-Amino-1-indanol, a useful synthetic intermediate that is currently being investigated for its potential use as a resolving agent for optical resolution. /// References: Saigo, K., et al. 3-Amino-1-indanol, Method for Synthesizing It From α -Phenyl- α -alanine, and Its Use As Resolving Agent for Optical Resolution. PCT Int. Appl. 2002020461. Mar 14, 2002

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------|---------------|---------------------|------|--------|--------|
| | (R,R)-trans-3-Amino-1-indanol | [403672-01-5] | | | | |