

(+)-6-Aminopenicillanic Acid-d3

Art. ID TRC-A618897-25MG
Unit 25 mg

Description

Application Notes: (+)-6-Aminopenicillanic Acid-d3 is the isotope labelled analog of (+)-6-Aminopenicillanic Acid. (+)-6-Aminopenicillanic Acid is the core structure in penicillins, obtained from the fermentation brew of the *Penicillium* mold. (+)-6-Aminopenicillanic Acid is used as the main starting block for the preparation of numerous semisynthetic penicillins. /// References: Zel'tser, I.Z. et al.: Gig. Truda Prof. Zabol., 5, 59 (1983)| Hamilton-Miller, J.M.T.: Int. J. Antimicrob. Agents, 31, 189 (2008)| Novikova, N.D. et al.: Mikol. Fitopatol., 7, 101 (1973)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(+)-6-Aminopenicillanic Acid-d3					