

### **(3S,4S)4-Amino-3-Pyrrolidinol Dihydrochloride**

Art. ID TRC-A629115-2.5G  
Unit 2.5 g

#### Description

Category: Heterocycles, Miscellaneous Reagents /// Application Notes: (3S,4S)4-Amino-3-Pyrrolidinol Dihydrochloride is used in the preparation of a novel alpha-mannosidase inhibitor. Also used in the synthesis of protected 4-hydroxypyrrolidines. /// References: Limberg, G. et al.: Synth., 1, 178 (1999)| Rodriguez-Rodriguez, J. et al.: Tetrahedron, 69, 5407 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(3S,4S)4-Amino-3-Pyrrolidinol Dihydrochloride	[220812-21-5]				