

1-Arachidoyl-sn-glycero-3-phosphocholine-d9

Art. ID TRC-A765097-25MG
Unit 25 mg

Description

Application Notes: 1-Arachidoyl-sn-glycero-3-phosphocholine-d9 is the isotope labelled analog of 1-Arachidoyl-sn-glycero-3-phosphocholine (A765095), a derivative of sn-Glycero-3-phosphocholine (G598700) and a biomarker in biological studies. /// References: Tang, C., et al.: J. Chromatogr. B, 879, 2095 (2011)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1-Arachidoyl-sn-glycero -3-phosphocholine-d9					