

## **ARP 100**

Art. ID                      TRC-A774058-100MG  
Unit                        100 mg

### Description

Application Note: ARP 100 is a potent, dual inhibitor of MMP-12/MMP-13 targets, which are a family of zinc dependent endopeptidases known to play key roles in extracellular matrix (ECM) breakdown disorders, such as the two main forms of arthritis, rheumatoid arthritis (RA) and osteoarthritis (OA).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	ARP 100	[704888-90-4]				