

Arsenazo-I Trisodium Hydrate

Art. ID TRC-A774170-2.5G
Unit 2.5 g

Description

Category: Amines, Aromatics, Pharmaceuticals, Intermediates and Fine Chemicals, Sulfur and Selenium Compounds /// Application Notes: Arsenazo-I is a reagent for the spectrophotometric determination of Al, Be, Ca, Th, U, V(IV), F- and SO₄²⁺, for Pd and Ru, for Ga and In, Arsenazo-I is a optical calcium indicator. /// References: Lamkin, E.G., et al.: Anal. Chem., 37, 1029 (1965)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenazo-I Trisodium Hydrate	[3547-38-4]				