

L-Aspartic Acid beta-(7-Amido-4-methylcoumarin)

Art. ID TRC-A788973-250MG

Unit 250 mg

Deliverydetails /not restricted

Description

Application Notes: L-Aspartic Acid beta-(7-Amido-4-methylcoumarin) is sensitive fluorogenic substrate used in lysosomal glycoasparaginase assays. This assay is good for diagnosis of aspartylglycosaaminuria.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Aspartic acid beta-(7 -amido-4-methylcoumarin)	[133628-73-6]				