

S-Acebutolol-d7

Art. ID TRC-A123812-100MG
Unit 100 mg

Description

Application Notes: S-Acebutolol-d7 is the isotope labelled analog of S-Acebutolol (A123812), the S-isomer of Acebutolol (A123800, HCl) which is a cardioselective beta-adrenergic blocker. Also antihypertensive, antianginal, and antiarrhythmic (class II). /// References: Cuthbert, O.A., et al.: Br. J. Pharmacol., 43, 639 (1971)| Singh, B.N., et al.: Drugs, 29, 531 (1985)| Foster, R.T., et al.: Anal. Profiles Drug Subs., 19, 1 (1990)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	S-Acebutolol-d7					