

## **N-Acetyl-D-histidine-3-15N**

Art. ID TRC-A178308-50MG  
Unit 50 mg

### Description

Application Notes: N-Acetyl-D-histidine-3-15N is the 15N labelled analogue of N-Acetyl-D-histidine (A178305) which was used in the preparation of oligonucleotide-ligand conjugates which may act as drug delivery and immunostimulant agents. /// References: Manoharan, M., et al.: Process for the preparation of oligonucleotide-ligand conjugates as drug delivery and immunostimulant agents, WO 2015006740, Jan 15, 2015|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-D-histidine-3- 15N					