

A-412997 Dihydrochloride

Art. ID TRC-A291708-250MG
Unit 250 mg
Deliverydetails UN2811 III E1

Description

Application Notes: The dopamine receptors play important role in cognition, memory, learning, and motor control (1). These receptors have been implicated as a therapeutic target for many psychiatric and neurological disorders. A-412997 Dihydrochloride is a selective agonist of the dopamine D4 receptor (2).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	A-412997 dihydrochlorid e	[1347744-96-0]				