

AGK 2

Art. ID TRC-A426600-250MG
Unit 250 mg
Deliverydetails /not restricted

Description

Category: Aromatics, Heterocycles, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: AGK 2 is a potent and selective histone deacetylase inhibitor or SIRT2 inhibitor, which may be used in the treatment of parkinsons disease, preventing cell death. /// References: Li, Y. et al.: Neurosci. Lett., 544, 36 (2013)| Rotili, D. et al.: J. Med. Chem., 55, 8193 (2012)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	AGK2	[304896-28-4]				