

DL-Alanine Methylamide

Art. ID TRC-A510590-2.5G
Unit 2.5 g

Description

Category: Amino Acids and Derivatives /// Application Notes: DL-Alanine Methylamide is a reagent used in the acetylation of DL-amino acid derivatives, which are potent new agents for the treatment of epilepsy. /// References: Conley, J.D., et al.: J. Med. Chem., 30, 567 (1987)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-Alanine Methylamide	[32012-12-7]				