

27-Alkyne Cholesterol

Art. ID TRC-A540300-250MG
Unit 250 mg

Description

Application Note: 27-Alkyne Cholesterol is an analog of alkyne cholesterol used for tracking cellular cholesterol metabolism and localization as it allows for manifold detection methods including mass spectrometry, thin-layer chromatog./fluorog., and fluorescence microscopy.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	27-Alkyne Cholesterol	[1527467-07-7]				