

(2-Amino-4-chlorophenyl)(tetrahydro-2H-pyran-3-yl)-methanone

Art. ID TRC-A604360-250MG

Unit 250 mg

Description

Application Notes: (2-Amino-4-chlorophenyl)(tetrahydro-2H-pyran-3-yl)-methanone is a reagent used as a building block in chemical synthesis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2-Amino-4-chlorophenyl) (tetrahydro-2H-pyran-3 -yl)-methanone	[1702437-90-8]				