

p-Aminomethyl Vorinostat

Art. ID TRC-A617070-25MG
Unit 25 mg

Description

Application Notes: p-Aminomethyl Vorinostat is a derivative of Vorinostat (S688700), a potent, selective, cell permeable histone deacetylase inhibitor (HDAC). /// References: Richon, V.M., et al.: Proc. Natl. Acad.Sci. USA, 95, 3003 (1988)| Munster, P.N., et al.: Cancer Res., 61, 8492 (2001)| Deroanne, C.F., et al.: Oncogene, 21, 247 (2002)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	p-Aminomethyl Vorinosta t	[1160823-16-4]				