

N-Arachidonoyltaurine-d4

Art. ID TRC-A765082-25MG
Unit 25 mg

Description

Category: Amines, Isotope Labelled Compounds, Sulfur and Selenium Compounds /// Application Notes:
N-Arachidonoyltaurine-d4 is the labeled analogue of N-Arachidonoyltaurine (A765080), a polyunsaturated lipid that inhibits the Cav3.3 current in T-type calcium channels. /// References: Cazade, M., et. al.: Mol. Pharmacol., 85, 218 (2014)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Arachidonoyltaurine-d					
	4					