

L-Argininamide-D-phenylalanyl-L-2-piperidinecarbonyl-N-(4-nitrophenyl)-dihydrochloride

Art. ID TRC-A769190-250MG
Unit 250 mg

Description

Application Note: L-Argininamide-D-phenylalanyl-L-2-piperidinecarbonyl-N-(4-nitrophenyl)-dihydrochloride has been used to study the mechanism by Which heparin proteoglycan modulates mast cell chymase activity and in the assay determination of heparin in blood.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Argininamide-D-phenyl alanyl-L-2-piperidineca rbonyl-N-(4-nitrophenyl)-dihydrochloride	[62354-65-8]				