

Azobenzene-4,4'-diboronic Acid

Art. ID TRC-A931005-100MG
Unit 100 mg

Description

Application Notes: Azobenzene-4,4'-diboronic Acid acts as a reagent for the synthesis and characterization of supramolecular gels based on boronate esters and imidazolyl donors. /// References: Luisier, N., et al.: Soft Matter, 12, 588 (2016)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Azobenzene-4,4'-diboronic Acid	[1902168-17-5]				