

(+)-cis,trans-Abscisic Acid-L-phenylalanine

Art. ID	TRC-A110050-25MG
Unit	25 mg
Deliverydetails	to be checked

Description

Application Note: (+)-cis,trans-Abscisic Acid-L-phenylalanine can be used in biological study for preparation and biological activity of amino acid conjugates of abscisic acid. Different sets of ABA-amino acids induced ABA-like responses in different plants. Different plants can also have hydrolyzing activity toward different ABA-amino acids and an ABA-amino acid may function as a species-selective pro-hormone of ABA.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(+)-cis,trans-Abscisic Acid-L-phenylalanine	[1285591-91-4]				