

N-Acetyl-beta-alanine-d3

Art. ID TRC-A168072-25MG
Unit 25 mg

Description

Application Notes: N-Acetyl-beta-alanine-d3 is the isotope labelled analog of N-Acetyl-beta-alanine (A168070), an abnormal amino acid metabolite formed in patients with isovaleric acidemia (IVA). /// References: Liu, Z., et al.: J. Med. Chem., 45, 631 (2002)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-beta-alanine-d					
	3					