

AHD-HN-13C3

Art. ID TRC-A434702-25MG
Unit 25 mg

Description

Category: Amines, Aromatics, Heterocycles, Isotope Labelled Compounds, Metabolites and Impurities ///

Application Notes: AHD-HN-13C3 is the isotope labelled analog of AHD-HN (A434700), a metabolite of 1-Amino Hydantoin (AHD, A611090) and 1-Formylnaphthalen-2-ol (F700825) which is a novel fluorescence agent. ///

References: Sheng, L., et al.: Food Addit. Contam., 30, 2114 (2013)| Du, N., et al.: J. Chromatogr. A, 1327, 90 (2014)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	AHD-HN-13C3					