

[[[4-(2-Amino-1,6-dihydro-6-oxo-9H-purin-9-yl)-2-butynyl]oxy]methyl]-phosphonic Acid

Art. ID TRC-A610080-100MG

Unit 100 mg

Description

Application Note: [[[4-(2-Amino-1,6-dihydro-6-oxo-9H-purin-9-yl)-2-butynyl]oxy]methyl]-phosphonic Acid is a rigid acyclonucleotide analog with antiviral activity.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	[[[4-(2-Amino-1,6-dihydro-6-oxo-9H-purin-9-yl)-2-butynyl]oxy]methyl]-phosphonic Acid	[159683-80-4]				