

(E)-2-Amino-4-(phosphonomethyl)hept-3-enoic Acid

Art. ID TRC-A627575-250MG
Unit 250 mg

Description

Application Notes: (E)-2-Amino-4-(phosphonomethyl)hept-3-enoic Acid is N-methyl-D-aspartate antagonist radioligand. It selectively binds to the NMDA receptor as an antagonist with high affinity and is the ligand of choice for labeling the NMDA receptor. /// References: Sills, M., et al.: Eur. J. Pharmacol., 192, 19 (1991)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(E)-2-Amino-4-(phosphonomethyl)hept-3-enoic Acid	[132472-31-2]				
	id					