

3-Aminopropyltriethoxy?silane-13C3

Art. ID TRC-A628787-25MG
Unit 25 mg

Description

Application Notes: 3-Aminopropyltriethoxy?silane-13C3 is the isotope labelled analog of 3-Aminopropyltriethoxy?silane (A628785), a reagent in the preparation of ultrastable organic dye nanoparticles in living cells. Also used in the preparation of monoclonal antibody-functionalized mesoporous silica nanoparticles packaged with Mifepristone (M343975), a progesterone receptor antagonist with partial agonist activity. Abortifacient. /// References: Xu, R. et al.: Biomat., 93, 38 (2016)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Aminopropyltriethoxy? silane-13C3					