

3-Amino-2,4,6-triiodobenzeneheptanoic Acid Ethyl Ester

Art. ID TRC-A630805-2.5G
Unit 2.5 g

Description

Category: Aliphatics, Glycerols, Fatty Acid Derivatives and Lipids, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: 3-Amino-2,4,6-triiodobenzeneheptanoic Acid Ethyl Ester is used to synthesize polyiodinated triglyceride analogs as potential computed tomography imaging agents for the liver. /// References: Weichert, J.P., et al.: J. Med. Chem., 38, 636 (1995)]

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Amino-2,4,6-triiodobenzeneheptanoic Acid Ethyl Ester	[161466-39-3]				