

Arachidonic Acid Methyl Ester-d3

Art. ID TRC-A765007-25MG
Unit 25 mg

Description

Category: Aliphatics, Glycerols, Fatty Acid Derivatives and Lipids, Isotope Labelled Compounds /// Application
Notes: Arachidonic Acid Methyl Ester-d3 is the isotope labelled analog of Arachidonic Acid Methyl Ester (A765005), the methyl ester derivative of Arachidonic Acid (A765000) that is shown to exhibit strong antimicrobial activity against oral microorganisms. Also used for preventing and treating gastric disorders. /// References: Huang, C.B., et al.: Arch. Oral Biol., 55, 555 (2010)| Choi, W., et al. PCT Int. Appl. WO 2010117194 A2 20101014. Oct 14, 2010

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arachidonic Acid Methyl Ester-d3					