

S-(+)-Arundic Acid Methyl Ester-d3

Art. ID TRC-A779527-10MG
Unit 10 mg

Description

Category: Isotope Labelled Compounds, Aliphatics, Chiral Reagents, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: S-(+)-Arundic Acid Methyl Ester-d3 is an isotope analog of S-(+)-Arundic Acid Methyl Ester. S-(+)-Arundic Acid Methyl Ester is a derivative of S-(+)-Arundic acid. S-(+)-Arundic acid (A779510) is an astrocyte modulating agent, in acute ischemic stroke. A therapeutic agent for the treatment of neurodegenerative diseases including Alzheimer's disease and Parkinson's disease. /// References: Uchida, Y., et al.: Brain Res., 481, 190 (1989), Tokunaga, M., et al.: Science, 277, 936 (1997), Sorbera, L., et al.: Drugs Future, 29, 441 (2004)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	S-(+)-Arundic Acid Methyl Ester-d3					