

L-Aspartyl-L-phenylalanine-d5

Art. ID TRC-A790012-25MG
Unit 25 mg

Description

Category: Amino Acids and Derivatives, Heterocycles, Impurities, Isotope Labelled Compounds, NeuroChemicals
Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: Labelled L-Aspartyl-L-phenylalanine
(A790010). L-Aspartyl-L-phenylalanine is an impurity of Thiamine (T344185). Labelled Thiamine Impurity B. ///
References: Aluwihare, L., et al.: Science, 308, 1007 (2005), Koch, B., et al.: Anal. Chem., 79, 1758 (2007),
Reemtsma, T., et al.: Environ. Sci. Technol., 42, 1430 (2008)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Aspartyl-L-phenylalanine-d5					