

S-(3-Acetylanisole-d3) Captopril

Art. ID TRC-A168532-100MG
Unit 100 mg

Description

Application Notes: S-(3-Acetylanisole-d3) Captopril is deuterium labeled S-(3-Acetylanisole) Captopril (A168530). S-(3-Acetylanisole) Captopril is a related compound of Captopril (C175750), which is a competitive inhibitor of angiotensin-converting enzyme (ACE) that converts angiotensin I (AT I) to angiotensin II (AT II). Captopril can be used to treat hypertension. // References: Rubin, B., et al.: Eur. J Pharmacol. 51, 377 (1978) | Ali, S. T. I., et al.: IRJP 6, 760 (2015)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	S-(3-Acetylanisole-d3)					
	Captopril					