

2-Acetylglutaric Acid Sodium Salt

Art. ID TRC-A178285-5G
Unit 5 g

Description

Category: Intermediates, Miscellaneous Reagents /// Application Notes: 2-Acetylglutaric Acid is used to synthesize coumarin beta-D-galactosides as novel fluorescent substrates for the detection of coliform bacteria in water. ///
References: Chilvers, K., et al.: J. Appl. Microbiol., 91, 1118 (2001)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Acetylglutaric Acid Sodium Salt					