

Adenosine Cyclic 2',3'-Phosphate 5'-Phosphate

Art. ID TRC-A280460-25MG
Unit 25 mg

Description

Application Notes: Adenosine Cyclic 2',3'-Phosphate 5'-Phosphate can be used in biological study for 2':3'-cyclic PAP, sulfuryl group acceptor in a reaction catalyzed by phenol sulfotransferase. It can also be used in biological study for effect of adenine nucleotides on myo-inositol-1,4,5-trisphosphate-induced calcium release. /// References: Lin, E. S., & Yang, Y. S.: Biochem. Bioph. Res. Co., 271, 818-822 (2000) | Missiaen, L., et al.: Biochem. J., 325, 661-666 (1997)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adenosine Cyclic 2',3'- Phosphate 5'-Phosphate	[4527-92-8]				