

D,L-Adrenochrome-d3

Art. ID TRC-A305252-25MG
Unit 25 mg

Description

Application Notes: D,L-Adrenochrome-d3 is the labelled analogue of D,L-Adrenochrome (A305250) which is the substance responsible for the red colours produced during the oxidation of adrenaline. /// References: Costa, V., et al.: Chem. Res. Toxicol., 20, 1183 (2007)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D,L-Adrenochrome-d3					