

N-Allylthiourea

Art. ID TRC-A572403-25G
Unit 25 g
Deliverydetails UN2811 III E1

Description

Application Notes: N-Allylthiourea is a reagent used in the synthesis of nano-composite materials and coatings used in the study of corrosion as well as in the preparation of carbon nanofiber-grafted highly porous polymeric beads for used as molecular filtration systems. /// References: Mishra, S. et al.: Chem. Eng. J., 313, 1142 (2017)| Balaji, J.. et al.: Prog. Org. Coat., 99, 463 (2016)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Allylthiourea	[109-57-9]				