

(2S)-2-Amino-3,3-di(methyl-d3)butanamide-4,4,4-d3

Art. ID TRC-A602692-500MG
Unit 500 mg

Description

Application Note: (2S)-2-Amino-3,3-di(methyl-d3)butanamide-4,4,4-d3 is used in the synthesis of deuterated analogs of the hepatitis C virus protease inhibitor boceprevir.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2S)-2-Amino-3,3-di(methyl-d3)butanamide-4,4,4-d3	[1092540-47-0]				