

3-(Aminoethyl)-5-methyl-4-hexenoic Acid

Art. ID TRC-A609025-100MG
Unit 100 mg

Description

Category: Miscellaneous Reagents /// Application Notes: 3-(Aminoethyl)-5-methyl-4-hexenoic Acid is used in the synthesis of aminoitransferase inhibitors. /// References: Lu, H. et al.: J. Med. Chem., 49, 7404 (2006)| Winter, D. et al.: J. Org. Chem., 75, 2610 (2010)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-(Aminoethyl)-5-methyl -4-hexenoic Acid					