

Amitifadine Lactone

Art. ID TRC-A633375-25MG
Unit 25 mg

Description

Application Note: Amitifadine Lactone is derived from 3,4-Dichlorophenylacetonitrile (D435435), which is an intermediate used in the synthesis of SR 48968, a potent, competitive and selective non-peptide antagonist of the neurokinin A (NK2) receptor¹ (1,2,3). It is also used to prepare tachykinin NK3 receptor antagonist SR 142801, (R)-N-[1-[3-[1-benzoyl-3-(3,4-dichlorophenyl)-3-piperidinyl]propyl]-4-phenyl-4-piperidinyl-N-methyl acetamide.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Amitifadine Lactone	[1643438-51-0]				