

(3R)-3-[4-Amino-3-(4-phenoxyphenyl)-1H-pyrazolo[3,4-d]pyrimidin-1-yl]-1-piperidinecarboxylic Acid Phenylmethyl Ester

Art. ID TRC-A634893-500MG

Unit 500 mg

Description

Application Note:

(3R)-3-[4-Amino-3-(4-phenoxyphenyl)-1H-pyrazolo[3,4-d]pyrimidin-1-yl]-1-piperidinecarboxylic Acid Phenylmethyl Ester is an intermediate in synthesizing Ibrutinib-d4 (I124971), a labelled analogue of Ibrutinib (I124970). Ibrutinib is a highly selective Bruton's tyrosine kinase (Btk) irreversible inhibitor.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(3R)-3-[4-Amino-3-(4-phenoxyphenyl)-1H-pyrazolo[3,4-d]pyrimidin-1-yl]-1-piperidinecarboxylic Acid Phenylmethyl Ester	[2012523-59-8]				