

## **(R)-Apomorphine Methanesulfinic Acid Salt**

Art. ID TRC-A727545-100MG  
Unit 100 mg

### Description

Application Notes: (R)-Apomorphine Methanesulfinic Acid Salt is a a dopamine (D1 and D2) receptor agonist derivative. May be used in the treatment of Parkinson's disease. /// References: Neumeyer, J.L., et al.: J. Pharm. Sci., 59, 1850 (1970), Ginos, J.Z., et al.: J. Med. Chem., 18, 1194 (1975), Muhtadi, F.J., et al.: Anal. Profiles Drug Subs., 20, 121 (1991)

| Text/Information | Analyte/Parameter                             | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------------------|------------|---------------------|------|--------|--------|
|                  | (R)-Apomorphine Methane<br>sulfinic Acid Salt |            |                     |      |        |        |