

(R)-Apremilast-d5

Art. ID TRC-A729707-100MG
Unit 100 mg

Description

Application Notes: Isotope labelled (R)-Apremilast is an enantiomer of Apremilast (A729700), an oral phosphodiesterase 4 inhibitor is used in the treatment of psoriatic arthritis. /// References: Gan, E. et al.: BioDrugs, 27, 359 (2013)| Pathan, E. et al.: Annal. Rheum. Dis., 72, 1475 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(R)-Apremilast-d5					