

AR-C133913XX-d7

Art. ID TRC-A766552-25MG
Unit 25 mg

Description

Application Notes: AR-C133913XX-d7 is the isotope labelled analog of AR-C133913XX (A766550), a metabolite of Ticagrelor (T437700) which is the first reversible oral P2Y₁₂ receptor antagonist that provides faster, greater, and more consistent ADP-receptor inhibition than Clopidogrel. Used in the treatment of acute coronary syndromes (ACS). /// References: Matetzky, S., et al.: Circulation, 109, 3171

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	AR-C133913XX-d7					