

(+)-cis,trans-Abscisic Acid-L-alanine

Art. ID TRC-A110065-25MG
Unit 25 mg
Deliverydetails to be checked

Description

Application Note: (+)-cis,trans-Abscisic Acid-L-alanine can be used in biological study for preparation and biological activity of amino acid conjugates of abscisic acid. Different sets of ABA-amino acids induced ABA-like responses in different plants. Different plants can also have hydrolyzing activity toward different ABA-amino acids and an ABA-amino acid may function as a species-selective pro-hormone of ABA.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(+)-cis,trans-Abscisic Acid-L-alanine	[1285591-80-1]				