

(S)-2-Acetoxypropionoyl Chloride-d3

Art. ID TRC-A167742-25MG
Unit 25 mg

Description

Application Note: (S)-2-Acetoxypropionoyl Chloride-d3 is an intermediate used in the synthesis of Iopamidol-d3 (I735602), a labeled non-ionic radiocontrast medium. Diagnostic aid (radiopaque medium).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(S)-2-Acetoxypropionoyl Chloride-d3					