

Adenosine-2',3'-cyclic Monophosphate Sodium Salt

Art. ID TRC-A280450-2.5MG
Unit 2.5 mg
Deliverydetails /not restricted

Description

Category: Heterocycles, Nucleotides, Bases and Related Reagents, Phosphorylating and Phosphitylating Agents ///
Appearance: White Solid /// Application Notes: Adenosine-2',3'-cyclic Monophosphate is an Adenosine (A280400)
derivative useful as anticancer agent. /// References: Touil, F., et al.: J. Pharm. Biomed. Anal., 40, 822 (2006)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adenosine 2':3'-cyclic monophosphate sodium sa lt	[37063-35-7]				