

Adenosine 5'-Triphosphate, 3'-phosphate, 1:1 mixture with Adenosine 5'-Triphosphate, 2'-phosphate

Art. ID TRC-A281815-25MG

Unit 25 mg

Description

Category: Carbohydrates and Derivatives, Heterocycles, Metabolites and Impurities, Nucleotides, Bases and Related Reagent, Phosphorylating and Phosphitylating Agents /// Application Notes: derivatives of Adenosine 5'-Triphosphate (A281810, Disodium Salt) which is used in cells as a coenzyme of intracellular energy transfer. It transports chemical energy within cells for metabolism. /// References: Knowles, J., et al.: Annu. Rev. Biochem., 49, 877 (1980)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adenosine 5'-Triphosphate, 3'-phosphate, 1:1 mixture with Adenosine 5'-Triphosphate, 2'-phosphate					