

(2S,4S,5S,7R)-Aliskiren Fumarate (2:1)

Art. ID TRC-A536020-50MG
Unit 50 mg

Description

Category: Amines, Chiral Reagents, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: (2S,4S,5S,7R)-Aliskiren Fumarate (2:1) is an isomer of Aliskiren Hemifumarate (A536001), an orally active, synthetic nonpeptide renin inhibitor. Antihypertensive. /// References: Lefevre, G., et al.: J. Immunoassay, 21, 65 (2000)| Allikments, K., et al.: Curr. Opin. Invest. Drugs, 3, 1479 (2002)| Wood| J.M., et al.: Biochem. Biophys. Res. Commun., 308, 698 (2003)| Grandman, A.H., et al.: Circulation, 111, 1012 (2005)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2S,4S,5S,7R)-Aliskiren	[1630036-82-6]				
	Fumarate (2:1)					