

(S,S)-trans-3-Amino-1-indanol

Art. ID TRC-A612110-250MG
Unit 250 mg

Description

Category: Amines, Aromatics /// Application Notes: (S,S)-trans-3-Amino-1-indanol is an isomer of 3-Amino-1-indanol, a useful synthetic intermediate that is currently being investigated for its potential use as a resolving agent for optical resolution. /// References: Saigo, K., et al. 3-Amino-1-indanol, Method for Synthesizing It From α -Phenyl- α -alanine, and Its Use As Resolving Agent for Optical Resolution. PCT Int. Appl. 2002020461. Mar 14, 2002

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(S,S)-trans-3-Amino-1-indanol	[403671-98-7]				