

(2R,3R,3aS,7aR,9R,10R,10aS,14aR)-2,9-bis(6-Amino-9H-purin-9-yl)dodecahydro-3,10-dihydroxydifuro[3,2-d:3',2'-j][1,3,7,9]tetraazacyclododecine-5,12(4H,6H)-dione

Art. ID TRC-A629830-25MG
Unit 25 mg
Deliverydetails to be checked

Description

Application Note:

(2R,3R,3aS,7aR,9R,10R,10aS,14aR)-2,9-bis(6-Amino-9H-purin-9-yl)dodecahydro-3,10-dihydroxydifuro[3,2-d:3',2'-j][1,3,7,9]tetraazacyclododecine-5,12(4H,6H)-dione, an amino acid derivative used in preparation of STING protein.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2R,3R,3aS,7aR,9R,10R,10aS,14aR)-2,9-bis(6-Amino-9H-purin-9-yl)dodecahydro-3,10-dihydroxydifuro[3,2-d:3',2'-j][1,3,7,9]tetraazacyclododecine-5,12(4H,6H)-dione	[2197188-04-6]				