

## **Arclofenin**

Art. ID                      TRC-A766585-1G  
Unit                            1 g

### Description

Application Notes: Arclofenin is derived from Nitrilotriacetic Acid (N490175), which is a chelating agent which forms coordination compound with metal ions. Nitrilotriacetic acid is used in complexometric titrations and as well as for protein isolation and purification in the His-tag method. /// References: Smolinska, B., et al: J. Chem. Tech. Biotech., 87, 1360 (2012)| Kaluza, L., et al.: Catalysis. Letters., 142, 969 (2012)| Tragardh, G., et al.: Desalination., 119, 21 (1998)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arclofenin	[87071-16-7]				